

# Initial S™: Single use sterile kits for hand surgery

## Initial S™ - Hand - Instrumentation kit content\*

### STERILE INSTRUMENTATION FOR Ø2.25 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC2.25H	Kit for Ø2.25 mm cannulated screws - Hand	Ø1.7 mm quick coupling drill bit - cannula Ø0.9 mm - L83 mm	1
		3 in 1 instrument for Ø2.25 mm screws	1
		Double guide Ø0.9 mm / Ø1.8 mm for Ø2.25 mm screws	1
		5.8 mm single use handle	1
		Pin Ø0.8 L150 mm	3

### STERILE INSTRUMENTATION FOR Ø2.6 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC2.6H	Kit for Ø2.6 mm cannulated screws - Hand	Ø2.0 mm quick coupling drill bit - cannula Ø1.1 mm - L83 mm	1
		3 in 1 instrument for Ø2.6 mm screws	1
		Double guide Ø1.1 mm / Ø2.1 mm for Ø2.6 mm screws	1
		5.8 mm single use handle	1
		Pin Ø1.0 L150 mm	3

### STERILE INSTRUMENTATION FOR Ø3.0 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC3.0H	Kit for Ø3.0 mm cannulated screws - Hand	Ø2.2 mm quick coupling drill bit - cannula Ø1.1 mm - L83 mm	1
		3 in 1 instrument for Ø3.0 mm screws	1
		Double guide Ø1.1 mm / Ø2.3 mm for Ø3.0 mm screws	1
		5.8 mm single use handle	1
		Pin Ø1.0 L150 mm	3

#### Indications :


The implants of the INITIAL S™ range (3.0 mm and smaller diameter) are intended for the fixation of bone fractures and osteotomies and for arthrodeses of foot and hand in adults.

#### Contra-indications :

- Serious vascular deterioration, bone devitalization,
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing a good fixation of the implants into the bone,
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

\* Sterile instrumentation for use with **Footmotion** screws range. Available separately.

For complete indications of the implants, please refer to the **Footmotion implants instructions for use**.



Example of instrumentation kit content for Ø3.0 mm screws.

# Initial S™ Implants references\*

\*Available separately in the Sterile Screw caddy.


## Ø2.25 MM CANNULATED SCREWS

Long thread	Short thread	Description
H0.9HFT2.25L10-ST	-	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 10 mm - STERILE
H0.9HFT2.25L12-ST	-	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 12 mm - STERILE
H0.9HFT2.25L14-ST	H0.9IFT2.25L14-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 14 mm - STERILE
H0.9HFT2.25L16-ST	H0.9IFT2.25L16-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 16 mm - STERILE
H0.9HFT2.25L18-ST	H0.9IFT2.25L18-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 18 mm - STERILE
H0.9HFT2.25L20-ST	H0.9IFT2.25L20-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 20 mm - STERILE
H0.9HFT2.25L22-ST	H0.9IFT2.25L22-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 22 mm - STERILE
H0.9HFT2.25L24-ST	H0.9IFT2.25L24-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 24 mm - STERILE
H0.9HFT2.25L26-ST	H0.9IFT2.25L26-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 26 mm - STERILE
H0.9HFT2.25L28-ST	H0.9IFT2.25L28-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 28 mm - STERILE
H0.9HFT2.25L30-ST	H0.9IFT2.25L30-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 30 mm - STERILE
H0.9HFT2.25L32-ST	H0.9IFT2.25L32-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 32 mm - STERILE
H0.9HFT2.25L34-ST	H0.9IFT2.25L34-ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 34 mm - STERILE

## Ø2.6 MM CANNULATED SCREWS

Long thread	Short thread	Description
H1.1HFT2.6L10-ST	-	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 10 mm - STERILE
H1.1HFT2.6L12-ST	-	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 12 mm - STERILE
H1.1HFT2.6L14-ST	H1.1IFT2.6L14-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 14 mm - STERILE
H1.1HFT2.6L16-ST	H1.1IFT2.6L16-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 16 mm - STERILE
H1.1HFT2.6L18-ST	H1.1IFT2.6L18-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 18 mm - STERILE
H1.1HFT2.6L20-ST	H1.1IFT2.6L20-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 20 mm - STERILE
H1.1HFT2.6L22-ST	H1.1IFT2.6L22-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 22 mm - STERILE
H1.1HFT2.6L24-ST	H1.1IFT2.6L24-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 24 mm - STERILE
H1.1HFT2.6L26-ST	H1.1IFT2.6L26-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 26 mm - STERILE
H1.1HFT2.6L28-ST	H1.1IFT2.6L28-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 28 mm - STERILE
H1.1HFT2.6L30-ST	H1.1IFT2.6L30-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 30 mm - STERILE
H1.1HFT2.6L32-ST	H1.1IFT2.6L32-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 32 mm - STERILE
H1.1HFT2.6L34-ST	H1.1IFT2.6L34-ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 34 mm - STERILE

## 2 TYPES OF CANNULATED SCREWS:



## Ø3.0 MM CANNULATED SCREWS

Long thread	Short thread	Description
H1.1HFT3.0L10-ST	-	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 10 mm - STERILE
H1.1HFT3.0L12-ST	-	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 12 mm - STERILE
H1.1HFT3.0L14-ST	H1.1IFT3.0L14-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 14 mm - STERILE
H1.1HFT3.0L16-ST	H1.1IFT3.0L16-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 16 mm - STERILE
H1.1HFT3.0L18-ST	H1.1IFT3.0L18-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 18 mm - STERILE
H1.1HFT3.0L20-ST	H1.1IFT3.0L20-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 20 mm - STERILE
H1.1HFT3.0L22-ST	H1.1IFT3.0L22-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 22 mm - STERILE
H1.1HFT3.0L24-ST	H1.1IFT3.0L24-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 24 mm - STERILE
H1.1HFT3.0L26-ST	H1.1IFT3.0L26-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 26 mm - STERILE
H1.1HFT3.0L28-ST	H1.1IFT3.0L28-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 28 mm - STERILE
H1.1HFT3.0L30-ST	H1.1IFT3.0L30-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 30 mm - STERILE
H1.1HFT3.0L32-ST	H1.1IFT3.0L32-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 32 mm - STERILE
H1.1HFT3.0L34-ST	H1.1IFT3.0L34-ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 34 mm - STERILE
H1.1HFT3.0L36-ST	-	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 36 mm - STERILE
H1.1HFT3.0L38-ST	-	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 38 mm - STERILE
H1.1HFT3.0L40-ST	-	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 40 mm - STERILE

The information presented in this brochure is intended to demonstrate a Newclip Technics product. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.



NEWCLIP TECHNICS  
PA de la Lande Saint Martin - 45 rue des Garotières  
44115 Haute Goulaine (France)  
Tél. : +33 (0)2 28 21 37 12 - Fax : +33 (0)2 40 63 68 37  
Mail : [orders@newcliptechnics.com](mailto:orders@newcliptechnics.com)  
[www.newcliptechnics.com](http://www.newcliptechnics.com)

NEWCLIP USA  
642 Larkfield Center  
Santa Rosa CA 95403, USA  
Phone : +1 707 230 5078  
Mail : [customerservice@newclipusa.com](mailto:customerservice@newclipusa.com)  
[www.newclipusa.com](http://www.newclipusa.com)

