

Initial S™: Single use sterile kits for foot surgery

Initial S™ - Foot - Instrumentation kit content*

STERILE INSTRUMENTATION FOR Ø2.25 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC2.25F	Kit for Ø2.25 mm cannulated screws - Foot	Ø1.7 mm quick coupling drill bit - cannula Ø0.9 mm - L75 mm	1
		3 in 1 instrument for Ø2.25 mm screws	1
		5.8 mm single use handle	1
		Pin Ø0.8 L80 mm	3

STERILE INSTRUMENTATION FOR Ø2.6 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC2.6F	Kit for Ø2.6 mm cannulated screws - Foot	Ø2.0 mm quick coupling drill bit - cannula Ø1.1 mm - L75 mm	1
		3 in 1 instrument for Ø2.6 mm screws	1
		5.8 mm single use handle	1
		Pin Ø1.0 L80 mm	3

STERILE INSTRUMENTATION FOR Ø3.0 MM CANNULATED SCREWS

Kit	Description	Content	Qty
KIT-SC3.0F	Kit for Ø3.0 mm cannulated screws - Foot	Ø2.2 mm quick coupling drill bit - cannula Ø1.1 mm - L80 mm	1
		3 in 1 instrument for Ø3.0 mm screws	1
		5.8 mm single use handle	1
		Pin Ø1.0 L80 mm	3

Indications :

The implants of the INITIAL S™ range are intended for the fixation of bone fractures and osteotomies and for arthrodeses of foot in adults.

Contra-indications :

- Serious vascular deterioration, bone devitalization,
- Pregnancy.
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing a good fixation of the implants into the bone,
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

STERILE INSTRUMENTATION FOR Ø2.0 MM SNAP OFF SCREWS

Kit	Description	content	Qty
KIT-S2.0	Kit for Ø2.0 mm Snap-off screws	Screwdriver for Ø2.0 mm Snap-off screws	1

* Sterile instrumentation for use with **Footmotion** screws range. Available separately.

Newclip Technics instruments are to be used only in trauma or orthopaedic surgery, when **Newclip Technics** implants for osteosynthesis are fitted.



Example instrument kit content for Ø3.0 mm screws.



NEWCLIP-TECHNICS




Initial S™ Implants references*

*Available separately in the **Sterile Screw caddy**.

Ø2.25 MM CANNULATED SCREWS

Ref.	Description
H0.9HFT2.25L10 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 10 mm - STERILE
H0.9HFT2.25L12 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 12 mm - STERILE
H0.9HFT2.25L14 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 14 mm - STERILE
H0.9HFT2.25L16 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 16 mm - STERILE
H0.9HFT2.25L18 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 18 mm - STERILE
H0.9HFT2.25L20 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 20 mm - STERILE
H0.9HFT2.25L22 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 22 mm - STERILE
H0.9HFT2.25L24 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 24 mm - STERILE
H0.9HFT2.25L26 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 26 mm - STERILE
H0.9HFT2.25L28 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 28 mm - STERILE
H0.9HFT2.25L30 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 30 mm - STERILE
H0.9HFT2.25L32 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 32 mm - STERILE
H0.9HFT2.25L34 - ST	Self-drilling self-compressive screw - Ø2.25 mm - cannulated Ø0.9 mm - L 34 mm - STERILE



Ø2.6 MM CANNULATED SCREWS

Ref.	Description
H1.1HFT2.6L10 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 10 mm - STERILE
H1.1HFT2.6L12 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 12 mm - STERILE
H1.1HFT2.6L14 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 14 mm - STERILE
H1.1HFT2.6L16 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 16 mm - STERILE
H1.1HFT2.6L18 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 18 mm - STERILE
H1.1HFT2.6L20 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 20 mm - STERILE
H1.1HFT2.6L22 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 22 mm - STERILE
H1.1HFT2.6L24 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 24 mm - STERILE
H1.1HFT2.6L26 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 26 mm - STERILE
H1.1HFT2.6L28 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 28 mm - STERILE
H1.1HFT2.6L30 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 30 mm - STERILE
H1.1HFT2.6L32 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 32 mm - STERILE
H1.1HFT2.6L34 - ST	Self-drilling self-compressive screw - Ø2.6 mm - cannulated Ø1.1 mm - L 34 mm - STERILE




Ø3.0 MM CANNULATED SCREWS

Ref.	Description
H1.1HFT3.0L10 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 10 mm - STERILE
H1.1HFT3.0L12 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 12 mm - STERILE
H1.1HFT3.0L14 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 14 mm - STERILE
H1.1HFT3.0L16 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 16 mm - STERILE
H1.1HFT3.0L18 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 18 mm - STERILE
H1.1HFT3.0L20 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 20 mm - STERILE
H1.1HFT3.0L22 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 22 mm - STERILE
H1.1HFT3.0L24 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 24 mm - STERILE
H1.1HFT3.0L26 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 26 mm - STERILE
H1.1HFT3.0L28 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 28 mm - STERILE
H1.1HFT3.0L30 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 30 mm - STERILE
H1.1HFT3.0L32 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 32 mm - STERILE
H1.1HFT3.0L34 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 34 mm - STERILE
H1.1HFT3.0L36 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 36 mm - STERILE
H1.1HFT3.0L38 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 38 mm - STERILE
H1.1HFT3.0L40 - ST	Self-drilling self-compressive screw - Ø3.0 mm - cannulated Ø1.1 mm - L 40 mm - STERILE



SNAP OFF SCREWS Ø2.0 MM

Ref.	Description
WST2.0L11 - ST	Snap-Off Weil screw - Ø2.0 mm - L 11 mm - STERILE
WST2.0L12 - ST	Snap-Off Weil screw - Ø2.0 mm - L 12 mm - STERILE
WST2.0L13 - ST	Snap-Off Weil screw - Ø2.0 mm - L 13 mm - STERILE
WST2.0L14 - ST	Snap-Off Weil screw - Ø2.0 mm - L 14 mm - STERILE
WST2.0L15 - ST	Snap-Off Weil screw - Ø2.0 mm - L 15 mm - STERILE



Sterile screw caddy



NEWCLIP TECHNICS
PA de la Lande Saint Martin - 45 rue des Garottières
44115 Haute Goulaine (France)
Tél. : +33 (0)2 28 21 37 12 - Fax : +33 (0)2 40 63 68 37
Mail : orders@newcliptechnics.com
www.newcliptechnics.com

NEWCLIP USA
642 Larkfield Center
Santa Rosa CA 95403, USA
Phone : +1 707 230 5078
Mail : customerservice@newclipusa.com
www.newclipusa.com

The information presented in this brochure is intended to demonstrate a Newclip Technics product. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.

